

Scope / Test® InstruSponge™ Usage Cross Reference Chart

VISN 10 Scope			Ruhof Test® InstruSponge™ Match		
SCOPE DESCRIPTION	MODEL NO.	CHANNEL DIA.	PROD NO.	SPONGE DIA.	BAG COLOR
Olympus Enteroscope	SIF100	2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-2T100	3.7/2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-P140	2.2mm	345TSPG2	2mm	YELLOW
Olympus Gastroscope	GIF-XQ140	2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-140	2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-Q140	2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-1T140	3.7mm	345TSPG4	4mm	PURPLE
Olympus Colonoscope	PCF-140L	3.2mm	345TSPG3	3mm	RED
Olympus Colonoscope	CF-140L	3.7mm	345TSPG4	4mm	PURPLE
Olympus Colonoscope	CF-Q140L	3.7mm	345TSPG4	4mm	PURPLE
Olympus Colonoscope	CF-1T140L	4.2mm	345TSPG4	4mm	PURPLE
Olympus Duodenoscope	JF-140	3.2mm	345TSPG3	3mm	RED
Olympus Duodenoscope	TJF-140	4.2mm	345TSPG4	4mm	PURPLE
Olympus Enteroscope	SIF-Q140	2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-XP160	2mm	345TSPG2	2mm	YELLOW
Olympus Gastroscope	GIF-160	2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-Q160	2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-2T160	3.7/2.8mm	345TSPG3	3mm	RED
Olympus Colonoscope	PCF-160AL	3.2mm	345TSPG3	3mm	RED
Olympus Colonoscope	CF-Q160AL/L	3.7mm	345TSPG4	4mm	PURPLE
Olympus Colonoscope	CF-Q160AL/AL	3.7mm	345TSPG4	4mm	PURPLE
Olympus Colonoscope	CF-2T160L	3.7/3.2mm	345TSPG3	3mm	RED
Olympus Duodenoscope	TJF-160F/VF	4.2mm	345TSPG4	4mm	PURPLE
Olympus Gastroscope	GIF-N180	2mm	345TSPG2	2mm	YELLOW
Olympus Gastroscope	GIF-Q180	2.8mm	345TSPG3	3mm	RED
Olympus Gastroscope	GIF-H180	2.8mm	345TSPG3	3mm	RED
Olympus Colonoscope	PCF-Q180AL	3.2mm	345TSPG3	3mm	RED
Olympus Colonoscope	CF-Q180AL	3.7mm	345TSPG4	4mm	PURPLE
Olympus Colonoscope	CF-H180AL	3.7mm	345TSPG4	4mm	PURPLE
Olympus Enteroscope	SIFQ180	2.8mm	345TSPG3	3mm	RED
Olympus Colonoscope	PCF H180AL/I	3.2mm	345TSPG3	3mm	RED
Olympus Radial Scope	GF-UE160-AL5	2.8mm	345TSPG3	3mm	RED
Olympus Convex Scope	GF-UC140P-AL5	2.8mm	345TSPG3	3mm	RED
Olympus Bronchoscope	BF-P180	2.0mm	345TSPG2	2mm	YELLOW
Olympus Bronchoscope	BF-Q180	2.0mm	345TSPG2	2mm	YELLOW
Olympus Bronchoscope	BF-1T180	3.0mm	345TSPG3	3mm	RED
Olympus Bronchoscope	BF-XT160	3.2mm	345TSPG3	3mm	RED
Olympus Duodenoscope	TJF-160	4.2mm	345TSPG4	4mm	PURPLE
Olympus Convex Scope	BF-UC160	2.0mm	345TSPG2	2mm	YELLOW
Olympus Gastroscope	GIF-N180	2.0mm	345TSPG2	2mm	YELLOW
Storz Adult Bronchoscope	11301BN1	2.3mm	345TSPG2	2mm	YELLOW
Storz Pediatric Bronchoscope	11302BD2	1.5mm	345TSPG1.5	1.5mm	BLACK
Storz Flexible Ureterscope	11278AU1	3.6fr	345TSPG1.5	1.5mm	BLACK
Storz Flexible Cystoscope	11272CU1	7fr	345TSPG2	2mm	YELLOW
Storz Flexible Cholescope	11292ADU1	3.6fr	345TSPG1.5	1.5mm	BLACK

• **Do not force Test® InstruSponge™ through the scope.** If there is any question regarding size, please call us and we will help you find the proper diameter Test® InstruSponge™.

• Dip the tip of the Test® InstruSponge™ in Sterile Water to help glide through scope channel easier and aid in the collection of contaminants.

• Please see Ruhof ATP Instruction Guide for complete Directions for Use.

